

HECTD2 antibody



Catalog Number: 111397

Product name

HECTD2 antibody

Immunogen

[Human HECTD2 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

HECTD2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 88 kDa

Preparation

This antibody was obtained by immunization of HECTD2 recombinant protein (Accession Number: NM_182765). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111397(HECTD2 Antibody) at dilution of 1:1000