

NT5M antibody

Catalog Number: 113397

Product name

NT5M antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NT5M Rabbit Polyclonal antibody. Positive IP detected in mouse skeletal muscle tissue. Positive WB detected in mouse brain tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 26-30 kDa

Preparation

This antibody was obtained by immunization of NT5M recombinant protein (Accession Number: NM_020201). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113397(NT5M antibody) at dilution of 1:600



IP Result of anti-NT5M (IP:Catalog No:113397, 4ug; Detection:Catalog No:113397 1:600) with mouse skeletal muscle tissue lysate 2200ug.