

# PHACS antibody



Catalog Number: 113797

## Product name

PHACS antibody

## Immunogen

[Human ACCS Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PHACS Rabbit Polyclonal antibody. Positive WB detected in MDA-MB-453s cells, HeLa cells, mouse testis tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 57 kDa

## Preparation

This antibody was obtained by immunization of PHACS recombinant protein (Accession Number: NM\_032592). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:113797(ACCS antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113797(ACCS antibody) at dilution of 1:50 (under 40x lens)