

Anti-HP1 alpha Rabbit antibody



Catalog Number: 168797

Product name

Anti-HP1 alpha Rabbit antibody

Dilutions

1:2000-1:10000

Specificity

1:50

Human, Mouse, Rat

1:50

Antibody description

1:20

HP1 alpha Rabbit monoclonal antibody

1:20


Validations

Preparation

Antigen: A synthetic peptide of human HP1 alpha

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Immunocytochemistry of HP1 alpha(green) in HeLa cells using HP1 alpha Rabbit mAb at dilution

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP



1/50, and DAPI(blue)

Western blot detection of HP1 alpha in 3T3, HeLa cell lysates using HP1 alpha Rabbit mAb(1:1000 diluted). Predicted band size: 22kDa. Observed band size: 22kDa.