

Human MAPKAPK5 (His & GST Tag) recombinant protein



Catalog Number: 500069

General Information

Gene Name Synonym

MAPK-activated protein kinase 5; MAPKAP kinase 5; MAPKAP-K5; MAPKAPK-5; MK-5; MK5; EC 2.7.11.1; p38-regulated/activated protein kinase; PRAK

Protein Construction

A DNA sequence encoding the human MAPKAPK5 isoform 1 (NP_003659.2) (Met 1-Gln 471) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

Kinase activity untested

Purity

> 76 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human MAPKAPK5/GST chimera consists of 708 amino acids and has a calculated molecular mass of 82 KDa. It migrates as an approximately 75KDa band in SDS-PAGE under reducing conditions.

Formulation

Supplied as sterile 20mM Tris, 500mM NaCl, pH 8.0, 10% glycerol

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human MAPKAPK5 Protein (His & GST Tag) SDS-PAGE