

# H7N1 HA (His Tag) recombinant protein



Catalog Number: 500088

## General Information

### Protein Construction

A DNA sequence encoding the Influenza A virus (A/turkey/Italy/4602/99(H7N1)) hemagglutinin (CAD38286.1) (Met1-Arg339), termed as HA1, was expressed with a polyhistidine tag at the C-terminus.

### Organism

H7N1

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per µg protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Asp 19

### Molecular Mass

The recombinant HA1 subunit of the Influenza A virus hemagglutinin

(A/turkey/Italy/4602/99(H7N1)) consists 332 amino acids and predicts a molecular mass of 36.5 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Influenza A H7N1 (A/turkey/Italy/4602/99)  
Hemagglutinin / HA1 Protein (His Tag)