

Human SRFBP1 (p49/STRAP) (His & GST Tag) recombinant protein



Catalog Number: 500090

General Information

Gene Name Synonym

SRF-dependent transcription regulation-associated protein; p49/STRAP

Protein Construction

A DNA sequence encoding the human SRFBP1(Q8NEF9) (Met1-Asp429) was expressed with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The secreted recombinant human SRFBP1

consists of 666 amino acids and predicts a molecular mass of 76.5 KDa. The apparent molecular mass of the protein is approximately 69 KDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, 10% gly, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human SRFBP1 Protein (His & GST Tag) SDS-PAGE