

# Human EDAR/Ectodysplasin A Receptor (His Tag) recombinant protein



Catalog Number: 500270

## General Information

### Gene Name Synonym

Anhidrotic ectodysplasin receptor 1; Downless homolog; EDA-A1 receptor; Ectodermal dysplasia receptor; Ectodysplasin-A receptor

### Protein Construction

A DNA sequence encoding the human EDAR (NP\_071731.1) extracellular domain (Met 1-Ile 189) was fused with a polyhistidine tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 96 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Glu 27

### Molecular Mass

The secreted recombinant human EDAR consists of 174 amino acids and has a predicted molecular mass of 19 kDa. The apparent molecular mass of the protein is approximately 30-40 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human EDAR / DL Protein (His Tag) SDS-PAGE