

H5N1 HA (His Tag) recombinant protein



Catalog Number: 500297

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/chicken/Yamaguchi/7/2004(H5N1)) hemagglutinin (BAD89305.1) (Met1-Gln530), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

Organism

H5N1

Expression Host

Human Cells

QC Testing

Purity

> (21.2+47.8+29.9) % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 17

Molecular Mass

The recombinant hemagglutinin of Influenza A

virus (A/chicken/Yamaguchi/7/2004(H5N1)) consists of 525 amino acids and predicts a molecular mass of 60 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H5N1 (A/chicken/Yamaguchi/7/2004) Hemagglutinin / HA Protein (His Tag)