

# Human CTRC/chymotrypsin C (His Tag) recombinant protein



Catalog Number: 500312

## General Information

### Gene Name Synonym

Caldecrin

### Protein Construction

A DNA sequence encoding the human CTRC (Q99895) (Met 1-Leu 268) was fused with a polyhistidine tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 97 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Cys 17

### Molecular Mass

The secreted recombinant human CTRC (pro

form) consists of 263 amino acids and has a predicted molecular mass of 29.3 kDa. It migrates as an approximately 36 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human CTRC / chymotrypsin C Protein (His Tag)  
SDS-PAGE