

# Cynomolgus TIM-3/HAVCR2 (Fc Tag) recombinant protein



Catalog Number: 500346

## General Information

### Protein Construction

A DNA sequence encoding the cynomolgus TIM3 (EHH54703.1) (Met1-Arg201) was expressed with the Fc region of human IgG1 at the C-terminus.

### Organism

Cynomolgus

### Expression Host

Human Cells

## QC Testing

### Purity

(70.9+26.4) % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Ser 22

### Molecular Mass

The recombinant cynomolgus TIM3 comprises 418 amino acids and has a calculated molecular mass of 46.4 KDa. The apparent molecular mass of the

protein is approximately 60.5 and 35.8 KDa in SDS-PAGE.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Cynomolgus TIM3 / HAVCR2 Protein (Fc Tag)  
SDS-PAGE