

## General Information

### Gene Name Synonym

Cytokine-responsive protein CR6; DNA damage-inducible transcript 2 protein

### Protein Construction

A DNA sequence encoding human GADD45G (O95257) (Met1-Glu159) was expressed.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Met

### Molecular Mass

The recombinant human GADD45G consists of 159

amino acids and predicts a molecular mass of 17.1 KDa. It migrates as an approximately 17 KDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human GADD45G / CR6 Protein SDS-PAGE