

General Information

Protein Construction

A DNA sequence encoding the N-terminal fragment (Met1-Arg362) of the influenza B hemagglutinin (influenza B virus (B/Brisbane/60/2008)) (ACN29383.1), termed as HA1, was fused with a polyhistidine tag at the C-terminus.

Organism

Influenza B

Expression Host

Human Cells

QC Testing

Purity

> 92 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 16

Molecular Mass

The secreted recombinant HA1 subunit of influenza A H5N1 HA (influenza B virus

(B/Brisbane/60/2008)) comprises 358 amino acids and has a predicted molecular mass of 39.1 kDa. As a result of glycosylation, it migrates as an approximately 66 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza B (B/Brisbane/60/2008) Hemagglutinin Protein (HA1 Subunit) (His Tag) SDS-PAGE