

Human SGK3/SGKL (His & GST Tag) recombinant protein



Catalog Number: 500432

General Information

Gene Name Synonym

Cytokine-independent survival kinase;
Serum/glucocorticoid-regulated kinase 3;
Serum/glucocorticoid-regulated kinase-like

Protein Construction

A DNA sequence encoding the human SGK3 (Q96BR1-1) (Met 1-Leu 496) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

No Kinase Activity

Purity

> 96 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human SGK3/GST chimera consists of 733 amino acids and has a calculated molecular mass of 85 kDa. It migrates as an approximately 68 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, 0.5mM PMSF, 10% gly, 1mM GSH, pH 7.4
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human SGK3 / SGKL Protein (His & GST Tag)
SDS-PAGE