

Mouse TNFSF12 (Fc Tag) recombinant protein



Catalog Number: 500517

General Information

Gene Name Synonym

TNF-related weak inducer of apoptosis; Tumor necrosis factor ligand superfamily member 12, membrane form; Tumor necrosis factor ligand superfamily member 12, secreted form

Protein Construction

A DNA sequence encoding the mouse Tnfsf12 (NP_035744.1) (Arg105-His249) was expressed with the Fc region of human IgG1 at the N-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> (83.8+3.6) % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu

Molecular Mass

The recombinant mouse Tnfsf12 consists 405 amino acids and predicts a molecular mass of 44.4 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Tnfsf12 protein SDS-PAGE