

Human IgG1 Fc recombinant protein



Catalog Number: 500583

General Information

Protein Construction

A DNA sequence encoding the human IgG1 Fc region (AAC82527.1) (Glu 99-Lys 330) (one aa mutation, 103 Cys/Ser) was expressed.

Organism

Human

Expression Host

CHO Cells

QC Testing

Purity

> 92 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 99

Molecular Mass

The recombinant human IgG1 Fc consists of 232 amino acids and has a predicted molecular mass of 26 kDa. As a result of glycosylation, the

apparent molecular mass of rhFc is approximately 34 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human IgG1 Fc Protein SDS-PAGE