

Human VEGF / VEGFA / VEGF165 recombinant protein



Catalog Number: 500706

General Information

Gene Name Synonym

Vascular permeability factor

Protein Construction

A DNA sequence encoding the human VEGF165 isoform (P15692-4) (Met1-Arg191) was expressed. Human and Cynomolgus VEGF165 sequences are identical.

Organism

Human

Expression Host

Human Cells

QC Testing

Activity

Measured in a cell proliferation assay using human umbilical vein endothelial cells (HUVEC). The ED₅₀ for this effect is typically 4-16 ng/mL.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 27

Molecular Mass

The recombinant human VEGF165 consists of 165 amino acids and predicts a molecular mass of 19.2 KDa. It migrates as an approximately 20 and 22 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 100 mM Glycine, 10 mM NaCl, pH 7.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human / Cynomolgus VEGF / VEGFA / VEGF165 Protein SDS-PAGE