

Canine HGF / Hepatocyte Growth Factor recombinant protein



Catalog Number: 500911

General Information

Gene Name Synonym

Hepatopoietin-A; Scatter factor; Hepatocyte growth factor alpha chain; Hepatocyte growth factor beta chain

Protein Construction

A DNA sequence encoding the canine HGF (NP_001002964.1) (Met1-Ser730) was expressed.

Organism

Canine

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 32

Molecular Mass

The recombinant canine HGF consists 699 amino acids and predicts a molecular mass of 79.9 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Canine HGF / Hepatocyte Growth Factor Protein
SDS-PAGE