

Mouse MBL-2 / MBL (His Tag) recombinant protein



Catalog Number: 501010

General Information

Gene Name Synonym

Mannan-binding protein; RA-reactive factor P28A subunit

Protein Construction

A DNA sequence encoding the mouse Mbl2 (NP_034906.1) (Glu19-Asp244) was expressed with a polyhistidine tag at the N-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant mouse Mbl2 consists 245 amino acids and predicts a molecular mass of 26.3 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse MBL-2 / MBL Protein (His Tag)