

Danio rerio (zebrafish) EFNB2A / Ephrin B2a (Fc Tag) recombinant protein



Catalog Number: 501068

General Information

Protein Construction

A DNA sequence encoding the zebrafish EFNB2A (O73874) (Met 1-Ala222) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Danio rerio (zebrafish)

Expression Host

Human Cells

QC Testing

Activity

Measured by its binding ability in a functional ELISA.

Immobilized human EphB4-His (Cat:502026) at 10 µg/ml (100 µl/well) can bind zebrafish EFNB2A-Fc, The EC₅₀ of zebrafish EFNB2A-Fc is 12.8-30.1 ng/ml.

Purity

(16.2+75.4) % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 25

Molecular Mass

The recombinant zebrafish EFNB2A consists of 439 amino acids and has a calculated molecular mass of 48.9 kDa. It migrates as an 97 and 58 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, PH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Danio rerio (zebrafish) EFNB2A / Ephrin B2a Protein (Fc Tag) SDS-PAGE