

# Mouse Chk1 (His & GST Tag) recombinant protein



Catalog Number: 501140

## General Information

### Gene Name Synonym

CHK1 checkpoint homolog; Checkpoint kinase-1

### Protein Construction

A DNA sequence encoding the mouse CHEK1 (O35280-1) (Met1-Thr476) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

### Organism

Mouse

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Activity

Kinase activity untested

### Purity

> 90 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu$ g of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Met

### Molecular Mass

The recombinant mouse CHEK1/GST chimera consists of 713 amino acids and has a calculated molecular mass of 82.2 kDa. The recombinant protein migrates as an approximately 78 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Supplied as sterile 20mM Tris, 500mM NaCl, pH 8.5, 3mM DTT, 10% gly

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse CHK1 / CHEK1 Protein (His & GST Tag)  
SDS-PAGE