

# Rat CD63 / Tspan-30 / Tetraspanin-30 (His Tag) recombinant protein



Catalog Number: 501149

## General Information

### Gene Name Synonym

Mast cell antigen AD1

### Protein Construction

A DNA sequence encoding the rat CD63 (P28648) (Ala103-Val203) was expressed with a polyhistidine tag at the N-terminus.

### Organism

Rat

### Expression Host

Human Cells

## QC Testing

### Purity

(25.6+25.0+35.9) % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

His

### Molecular Mass

The recombinant rat CD63 comprises 121 amino

acids and has a predicted molecular mass of 13.8 kDa. The apparent molecular mass of the protein is approximately 24, 19, 15 kDa in SDS-PAGE under reducing conditions due to glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Rat CD63 / Tspan-30 / Tetraspanin-30 Protein (His Tag) SDS-PAGE