

Canine ANGPTL7/Angiopoietin-like 7 (His Tag) recombinant protein



Catalog Number: 501357

General Information

Protein Construction

A DNA sequence encoding the canine ANGPTL7 (XP_852356.1) (Met1-Pro344) was expressed with a C-terminal polyhistidine tag.

Organism

Canine

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 27

Molecular Mass

The recombinant canine ANGPTL7 comprises 329 amino acids and has a predicted molecular mass of 38.2 kDa. The apparent molecular mass of the

protein is approximately 37 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Canine ANGPTL7 / Angiopoietin-like 7 Protein (His Tag) SDS-PAGE