

# Mouse TGF-beta 2 (His Tag) recombinant protein



Catalog Number: 501424

## General Information

### Gene Name Synonym

Latency-associated peptide

### Protein Construction

A DNA sequence encoding the mouse Tgfb2 (NP\_033393.2) (Met1-Ser414) was expressed with a polyhistidine tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> (34.0+44.9+12.0)% as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per µg protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Ser 20

### Molecular Mass

The recombinant mouse Tgfb2 consists 406 amino acids and predicts a molecular mass of 46.9 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse TGF-beta 2 / TGFB2 Protein (His Tag)