

Rat IFN-alpha / IFNA1 / IFN (His Tag) recombinant protein



Catalog Number: 501742

General Information

Protein Construction

A DNA sequence encoding the rat IFNA1 (P05011) (Met1-Ser192) was expressed, fused with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Activity

Measured in antiviral assays using L929 cells infected with vesicular stomatitisvirus (VSV). The ED₅₀ for this effect is typically 4-20 ng/mL.

Purity

> 99 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Cys 24

Molecular Mass

The recombinant rat IFNA1 comprises 180 amino acids and predicts a molecular mass of 20.9 kDa. The apparent molecular mass of the rat IFNA1 is approximately 20-25 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat IFN-alpha / IFNA1 / IFN Protein (His Tag)
SDS-PAGE