

Human IL18R1 (His & Fc Tag) recombinant protein



Catalog Number: 501773

General Information

Gene Name Synonym

CD218 antigen-like family member A; CDw218a; IL1 receptor-related protein

Protein Construction

A DNA sequence encoding the extracellular domain of human IL18R1 (NP_003846.1) (Met 1-Arg 329) was fused with the C-terminal polyhistidine-tagged Fc region of human IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to bind cynomolgus IL-18 in a functional ELISA.

Purity

> 70 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 19

Molecular Mass

The recombinant human IL18R1/Fc is a disulfide-linked homodimer. The reduced monomer comprises 559 amino acids and has a predicted molecular mass of 64 kDa. As a result of glycosylation, the apparent molecular mass of rh IL18R1/Fc monomer is approximately 95-110 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human IL18R1 / CD218a Protein (His & Fc Tag)
SDS-PAGE