

Human NPY/Pro-neuropeptide Y (Fc Tag) recombinant protein



Catalog Number: 501950

General Information

Gene Name Synonym

cDNA, FLJ92442, Homo sapiens neuropeptide Y (NPY)

Protein Construction

A DNA sequence encoding the human NPY (NP_000896.1) (Tyr29-Tyr64) was expressed with the Fc region of human IgG1 at the N-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu

Molecular Mass

The recombinant human NPY consists 296 amino acids and predicts a molecular mass of 32.7 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human NPY / Neuropeptide Y Protein (Fc Tag)
SDS-PAGE