

Cynomolgus TNFR2/CD120b/TNFRSF1B (His Tag) recombinant protein



Catalog Number: 502021

General Information

Gene Name Synonym

Uncharacterized protein

Protein Construction

A DNA sequence encoding the cynomolgus TNFRSF1B (F7EAE4) (Met1-Asp257) was expressed, fused with a polyhistidine tag at the C-terminus.

Organism

Cynomolgus

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to inhibit the TNF α mediated cytotoxicity in the L-929 mouse fibroblast cells in the presence of the metabolic inhibitor actinomycin D.
The ED₅₀ for this effect is typically 0.15-0.45 μ g/mL in the presence of 1 ng/mL of cynoTNF α .

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 23

Molecular Mass

The recombinant cynomolgus TNFRSF1B comprises 246 amino acids and has a calculated molecular mass of 26.5 KDa. The apparent molecular mass of it is approximately 38 KDa respectively in SDS-PAGE.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Cynomolgus TNFR2 / CD120b / TNFRSF1B
Protein (His Tag) SDS-PAGE