

Human SNAP-alpha (His Tag) recombinant protein



Catalog Number: 502379

General Information

Protein Construction

A DNA sequence encoding the human NAPA (NP_003818.2) (Met1-Arg295) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human NAPA consists 310 amino acids and predicts a molecular mass of 35 kDa.

Formulation

Lyophilized from sterile 50 mM Tris, 10 % Glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human NAPA / alpha-SNAP Protein (His Tag)
SDS-PAGE