

Rat IL7 / interleukin-7 recombinant protein



Catalog Number: 502444

General Information

Protein Construction

A DNA sequence encoding the rat IL7 (Q91Y32) (Asp26-Ile154) was expressed and purified with an initial Met.

Organism

Rat

Expression Host

E. coli

QC Testing

Activity

Measured by its ability to bind rat IL7R-mFc (Cat:501881) in a functional ELISA.

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant rat IL7 consists of 130 amino acids and has a predicted molecular mass of 15 kDa. The apparent molecular mass of the recombinant protein is approximately 17 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat IL7 / interleukin 7 Protein SDS-PAGE