

# Human p38 alpha/MAPK14 (His Tag) recombinant protein



Catalog Number: 502461

## General Information

### Protein Construction

A DNA sequence encoding the full length of human MAPK14 isoform 2 (NP\_620581.1) (Met 1-Ser 360) was fused with a polyhistidine tag at the carboxy-terminus.

### Organism

Human

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Activity

No Kinase Activity

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu$ g of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Met 1

### Molecular Mass

The recombinant human MAPK14 consists of 370 amino acids and predicts amolecular mass of 42.8 kDa which is also estimated by SDS-PAGE.

### Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 7.4, 10% glycerol

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human p38 alpha / MAPK14 Protein (His Tag)  
SDS-PAGE