

Human AIM2 (GST Tag) recombinant protein



Catalog Number: 502867

General Information

Gene Name Synonym

Absent in melanoma 2

Protein Construction

A DNA sequence encoding the human AIM2 (NP_004824.1) (Met 1-Thr 343) was fused with the GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

Measured by its ability to inhibit the proliferation of MCF7 human breast adenocarcinoma cells [Karey, K.P. et al. (1988) Cancer Research 48:4083.]. The ED₅₀ for this effect is typically 4-7 µg/ml.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human AIM2/GST chimera consists of 567 amino acids and predicts a molecular mass of 65.2 kDa as estimated in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 50mM Tris, 1M NaCl, 0.5mM PMSF, 5mM GSH, pH 8.0
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human AIM2 Protein (GST Tag) SDS-PAGE