

# Rat NPC2 (His Tag) recombinant protein



Catalog Number: 503090

## General Information

### Gene Name Synonym

Niemann-Pick disease, type C2; Protein Npc2

### Protein Construction

A DNA sequence encoding the rat Npc2 (CAD56199.1) (Met1-Gly152) was expressed with a polyhistidine tag at the C-terminus.

### Organism

Rat

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Glu 23

### Molecular Mass

The recombinant rat Npc2 consists 141 amino acids and predicts a molecular mass of 15.9 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Rat Niemann-Pick disease type C2 / NPC2 Protein (His Tag)