

Human ALOX15B / 15 Lipoxygenase 2 (His & GST Tag) recombinant protein



Catalog Number: 503177

General Information

Gene Name Synonym

15-lipoxygenase 2; Arachidonate 15-lipoxygenase type II; Linoleate 13-lipoxygenase 15-LOb

Protein Construction

A DNA sequence encoding the human ALOX15B (AAH35217.1)(Met1-Ile676) was expressed the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The secreted recombinant human ALOX15B

consists of 913 amino acids and predicts a molecular mass of 103.7 KDa. The apparent molecular mass of the protein is approximately 93 KDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, 10% glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human ALOX15B / 15 Lipoxygenase 2 Protein (His & GST Tag) SDS-PAGE