

Human IL-8/CXCL8 (His Tag) recombinant protein



Catalog Number: 503360

General Information

Gene Name Synonym:

Protein Construction:

A DNA sequence encoding the human CXCL8 (NP_000575.1) (Ser28-Ser99) was expressed with a polyhistidine tag at the C-terminus.

Source:Human

Expression Host:Yeast

QC Testing

Activity:

Purity:

> 95 % as determined by SDS-PAGE.

Endotoxin:

Please contact us for more information.

Stability:

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal:

Ser 28

Molecular Mass:

The recombinant human CXCL8 consists of 82

amino acids and predicts a molecular mass of 9.8 kDa.

Formulation:

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage:

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution:

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE:



Human IL-8 / CXCL8 Protein (His Tag)