

Human IFNGR1/CD119 (His Tag) recombinant protein



Catalog Number: 503455

General Information

Gene Name Synonym

CDw119

Protein Construction

A DNA sequence encoding the extracellular domain (Met 1-Gly 245) of human IFN- γ R1 pre-protein (NP_000407.1) was expressed, fused with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to inhibit rhIFN γ mediated protection of WISH cells infected with vesicular stomatitis virus(VSV).

The ED₅₀ for this effect is typically 2-8 μ g/mL.

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 18

Molecular Mass

The recombinant human IFN- γ R1 consists of 239 amino acids after removal of the signal peptide and predicts a molecular mass of 27.3 kDa. By SDS-PAGE under reducing conditions, the apparent molecular mass of rhIFN- γ R1 is approximately 40-35 kDa due to the glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human IFNGR1 / CD119 Protein (His Tag) SDS-PAGE