

# Mouse PGF (Fc Tag) recombinant protein



Catalog Number: 503524

## General Information

### Protein Construction

A DNA sequence encoding the mouse PGF (P49764) (Met 1-Pro 158) was fused with the Fc region of human IgG1 at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> 76 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu$ g of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Val 19

### Molecular Mass

The recombinant mouse PGF/Fc is a disulfide-linked homodimer. The reduced monomer comprises 381 amino acids and has a calculated molecular mass of 42.9 kDa. As a result of

glycosylation, the apparent molecular mass of the monomer is approximately 52 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse PIGF / PLGF Protein (Fc Tag) SDS-PAGE