

Mouse CD32/Fc gamma RII/FCGR2 (His & AVI Tag) recombinant protein



Catalog Number: 503525

General Information

Protein Construction

A DNA sequence encoding the extracellular domain (Met 1-Arg 217) of mouse FGGR2B (NP_001070657.1) as fused with a c-terminal polyhistidine tagged AVI tag at the C-terminus. The expressed protein was biotinylated *in vivo* by the Biotin-Protein ligase (BirA enzyme) which is co-expressed.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Activity

1. Measured by its ability to bind human IgG1 in a functional ELISA.
2. Labeling ratio of biotin to protein: 0.5

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Thr 40

Molecular Mass

The secreted recombinant mouse FGGR2B consists of 212 amino acids and has a predicted molecular mass of 23.9 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rm FGGR2B is approximately 38-43 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CD32 / FCGR2B Protein (His & AVI Tag), Biotinylated SDS-PAGE