

# Human megakaryocyte potentiating factor (Fc Tag) recombinant protein



Catalog Number: 503754

## General Information

### Gene Name Synonym

CAK1 antigen; Pre-pro-megakaryocyte-potentiating factor; Megakaryocyte-potentiating factor; Mesothelin, cleaved form

### Protein Construction

A DNA sequence encoding the human MSLN (Q13421-2) (Met1-Arg286) was expressed, fused with the Fc region of human IgG1 at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 99 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Leu 37

### Molecular Mass

The recombinant human MSLN/Fc comprises 491 amino acids and has a predicted molecular mass of 53.7 kDa. The apparent molecular mass of the protein is approximately 56 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human MSLN / Mesothelin Protein (Fc Tag) SDS-PAGE