

Mouse CaM Kinase 4/CaMKIV recombinant protein



Catalog Number: 503827

General Information

Gene Name Synonym

CaM kinase-GR

Protein Construction

A DNA sequence encoding the mouse CAMK4(P08414) (Met1-Tyr469) was expressed and purified with two additional amino acids (Gly & Pro) at the N-terminus.

Organism

Mouse

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

Kinase activity untested

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly

Molecular Mass

The recombinant mouse CAMK4 consists of 471 amino acids and predicts a molecular mass of 52.7 KDa. It migrates as an approximately 55 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, 10% glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CAMK4 / CaMKIV Protein SDS-PAGE