

General Information

Gene Name Synonym

Uncharacterized protein

Protein Construction

A DNA sequence encoding the cynomolgus DPP4 (F6VRB0) (Asp34-Pro766) was expressed with two additional amino acids (Gly & Pro) at the N-terminus.

Organism

Cynomolgus

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly

Molecular Mass

The recombinant cynomolgus DPP4 is a disulfide-

linked homodimer. The reduced monomer comprises 735 amino acids and has a calculated molecular mass of 84.6 KDa. The apparent molecular mass of the protein is approximately 92-102 KDa in SDS-PAGE.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Cynomolgus DPP4 / DPPIV / CD26 Protein SDS-PAGE