

Mouse GPNMB/Osteoactivin (His Tag) recombinant protein



Catalog Number: 503904

General Information

Gene Name Synonym

DC-HIL; Dendritic cell-associated transmembrane protein; Osteoactivin

Protein Construction

A DNA sequence encoding the mouse GPNMB (Q99P91) (Met1-Asn502) was expressed with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Lys 23

Molecular Mass

The recombinant mouse GPNMB comprises 490 amino acids and has a predicted molecular mass of 54.9 kDa. The apparent molecular mass of the protein is approximately 96 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse GPNMB / Osteoactivin Protein (His Tag)
SDS-PAGE