

Mouse JAML/AMICA1 (Fc Tag) recombinant protein



Catalog Number: 504025

General Information

Gene Name Synonym

Dendritic cell-specific protein CREA7

Protein Construction

A DNA sequence encoding the mouse (NP_001005421.3) (Met1-Leu281) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> (7.3+91.7) % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 23

Molecular Mass

The recombinant mouse consists 500 amino acids and predicts a molecular mass of 56.5 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse JAML / AMICA1 Protein (Fc Tag) SDS-PAGE