

Human ILT6 / LILRA3 (His Tag) recombinant protein



Catalog Number: 504071

General Information

Gene Name Synonym

CD85 antigen-like family member E;
Immunoglobulin-like transcript 6; Leukocyte
immunoglobulin-like receptor 4; Monocyte
inhibitory receptor HM43/HM31

Protein Construction

A DNA sequence encoding the human LILRA3
(AAH28208.1) (Met1-Glu439) was expressed with
a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by
the LAL method

Stability

Samples are stable for up to twelve months from
date of receipt at -70°C

Predicted N terminal

Gly 24

Molecular Mass

The recombinant human LILRA3 consists of 427
amino acids and predicts a molecular mass of 46.5
KDa. It migrates as an approximately 57-66 KDa
band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as
protectants before lyophilization.
2. Please contact us for any concerns or special
requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C .
It is recommended that the protein be aliquoted
for optimal storage. Avoid repeated freeze-thaw
cycles.

Reconstitution

Adding sterile water, prepare a stock solution of
0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human ILT6 / LILRA3 Protein (His Tag) SDS-
PAGE