

## General Information

### Gene Name Synonym

Apolipoprotein A1; Proapolipoprotein A-I;  
Truncated apolipoprotein A-I; Apolipoprotein A-I(1-242)

### Protein Construction

A DNA sequence encoding the mature form of human APOA1 (P02647) (Asp25-Gln267) was expressed with a polyhistidine tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

His

### Molecular Mass

The recombinant human APOA1 consists of 263 amino acids and predicts a molecular mass of 30.7 KDa. It migrates as an approximately 27-31 KDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, 10% Glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human APOA1 / ApoAI Protein (His Tag) SDS-PAGE