

Mouse CD155/PVR (His & Fc Tag) recombinant protein



Catalog Number: 504218

General Information

Gene Name Synonym

Poliovirus receptor; Tage4 receptor

Protein Construction

A DNA sequence encoding the mouse CD155 (NP_081790.1) precursor (Met 1-Arg 345) was fused with the C-terminal polyhistidine-tagged Fc region of human IgG1 at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Activity

Measured by its binding ability in a functional ELISA. Immobilized mouse DNAM1 at 2 µg/ml (100 µl/well) can bind mouse CD155 with a linear ranger of 6.4-160 ng/ml.

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly 28

Molecular Mass

The secreted recombinant mouse CD155/Fc is a disulfide-linked homodimer after removal of the signal peptide. The reduced monomer comprises 566 amino acids with a predicted molecular mass of 63 kDa. As a result of glycosylation, it migrates as an approximately 95-105 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CD155 / PVR Protein (His & Fc Tag) SDS-PAGE