

# H10N8 HA (His Tag) recombinant protein



Catalog Number: 504653

## General Information

### Protein Construction

A DNA sequence encoding the Influenza A virus (A/Jiangxi-Donghu/346/2013(H10N8)) hemagglutinin (EPI497477) (Met1-Asp524) was expressed with a C-terminal polyhistidine tag.

### Organism

H10N8

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Activity

Measured by its ability to agglutinate guinea pig red blood cells.  
HA titer is 0.04-0.2 µg/mL for 1% GRBC.

### Purity

> 90 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Leu 17

### Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/Jiangxi-Donghu/346/2013(H10N8)) comprises 519 amino acids and has a predicted molecular mass of 57.9 kDa.

### Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



HA protein SDS-PAGE