

Rat CEACAM1/CD66a (His Tag) recombinant protein



Catalog Number: 504725

General Information

Protein Construction

A DNA sequence encoding the rat CEACAM1 (Met1-Ser422) was expressed with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 35

Molecular Mass

The recombinant rat CEACAM1 comprises 399 amino acids and predicts a molecular mass of 44.9 kDa. The apparent molecular mass of the recombinant protein is approximately 66-76 kDa

in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat CEACAM1 / CD66a Protein (His Tag) SDS-PAGE