

Human AOC1/ABP1/DAO1 (His tag) recombinant protein



Catalog Number: 518204

General Information

Gene Name Synonym

Amiloride-sensitive amine oxidase [copper-containing] DAO Diamine oxidase Amiloride-binding protein 1 Amine oxidase copper domain-containing protein 1 Histaminase Kidney amine oxidase KAO AOC1 ABP1, DAO1

Protein Construction

A DNA sequence encoding the Human AOC1/ABP1/DAO1 Glu20-Val751 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human AOC1/ABP1/DAO1 (His tag) recombinant protein