

# Human LOXL2 (Avi & His tag) recombinant protein



Catalog Number: 518209

## General Information

### Gene Name Synonym

Lysyl oxidase homolog 2; Lysyl oxidase-like protein 2; Lysyl oxidase-related protein 2; Lysyl oxidase-related protein WS9-14; LOXL2

### Protein Construction

A DNA sequence encoding the Human LOXL2 Gln26-Gln774 is expressed fused with a Avi-His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 50mM Sodium borate, 10mM CaCl<sub>2</sub>, 1.2M Urea, pH8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human LOXL2 (Avi & His tag) recombinant protein